

Activation and Startup Time

When the flow meter is activated, the startup process will begin, indicated by a white screen. After about 20 seconds, the loading animation will begin, followed by the company logo. Startup time typically takes 25-30 seconds after activation, at which time the default measurement screen will appear (see **Fig. 1** below)

Default Screen

The default screen of the flow meter is divided into four distinct display areas: flow, temperature, atmospheric pressure, and back pressure. Each of these is separated into their own colored box for easy distinction. At the top right of the screen is a gear icon, which, when pressed, will bring the user into the settings screen (Fig. 2) .




Fig. 1
Default screen




Fig. 2
Settings Screen

Changing Settings

There are two main methods for changing a setting of any given measurement. The first method is broken down into these steps:

1. Tap the gear icon in the top right of the default screen to enter the settings menu
2. Tap the corresponding name or icon of the setting you wish to change
3. Select any necessary settings in the appropriate submenu
4. Tap "Save" or "Cancel" to return to the settings menu
5. Tap "Back" or the arrow icon on the top right of the settings menu to return to the default screen

Mass Flow Rate



Tap once for settings menu

70.34 g/min
Air

Back Pressure 8699.1 mm Hg

Atm. Pressure 7353.8 millibar abs

Temperature 68.2 °C

 Device Info	<h1>Settings</h1>	Back 
Flow 	Back Pressure 	Atmospheric Pressure 
Temperature 	Response Time 	Set Zero 0.0 

Settings Screen

Settings

Medium	Units	
<input type="text" value="Air"/> Oxygen Air+Oxygen	<input type="text" value="scfh"/> lbm/min kg/min g/sec g/min mg/sec	
<input type="text" value="Oxygen %"/>		
<input type="button" value="Save"/>	<input type="button" value="Cancel"/>	<input type="button" value="Main Display"/>

Flow Settings

Alternatively, the settings for a single measurement can be adjusted as follows:

1. For 2 seconds, tap and hold down anywhere in the box of the measurement you wish to change
2. When the settings screen for that measurement appears, make the desired changes to the setting
3. Tap "Save" or "Cancel" to return to the default screen

Mass Flow Rate 

70.34 g/min
Air

8699.1 mm Hg

7353.8 millibar abs

68.2 °C

 Tap and hold for two seconds

Settings

Medium	Units
<div style="border: 1px solid black; padding: 2px;">Air Oxygen Air+Oxygen</div>	<div style="border: 1px solid black; padding: 2px;">scfh lbm/min kg/min g/sec g/min mg/sec</div>
<div style="border: 1px solid gray; padding: 2px; width: fit-content; margin: 10px auto;">Oxygen %</div>	
<div style="display: flex; justify-content: space-around; margin-top: 10px;"><div style="background-color: #0056b3; color: white; padding: 10px 20px; border-radius: 5px;">Save</div><div style="background-color: #0056b3; color: white; padding: 10px 20px; border-radius: 5px;">Cancel</div><div style="background-color: #4caf50; color: white; padding: 10px 20px; border-radius: 5px;">Main Display</div></div>	

The settings for Flow appear

Mass Flow Rate

70.34

g/min
Air

Back Pressure	8699.1 mm Hg
Atm. Pressure	7353.8 millibar abs
Temperature	68.2 °C

Tap and hold
for two seconds

Settings

Back Pressure

kPa	▲
psi	▼
bar	
millibar	
in Hg	▼

Save	Cancel	Main Display
------	--------	--------------


Back Pressure setting appears

Main Display setting

One setting that all four measurements share is the "Main Display" setting, located in the bottom right of each setting screen next to the "Cancel" button. When highlighted green, this indicates that the selected measurement will be at the top of the default display, taking up the most amount of space (it will become the "main display"). An example of changing between using Flow as the main display and Atmospheric Pressure is seen below.

Mass Flow Rate	
70.34	g/min Air
Back Pressure	8699.1 mm Hg
Atm. Pressure	7353.8 millibar abs
Temperature	68.2 °C

Flow selected as main display



Atmospheric Pressure selected as main display

Air + Oxygen specific setting

There is a specific setting when the "Air + Oxygen" medium is selected that is only selectable with this medium called "Oxygen %" -- this option will be grayed out when any other medium is selected.

This setting specifies what percentage of oxygen makes up the specific air + oxygen mixture. Valid oxygen percentages are 22% to 99% (21% or 100% would simply require selecting the "Air" or "Oxygen" mediums, specifically) and is selected by tapping the corresponding +1, +10, -1, or -10 buttons. The +10 and -10 buttons will not work if the resulting value is higher than 99% or lower than 1%, and so the +1 or -1 buttons must be used in these cases.

Settings

Medium

Units

Air
Oxygen
Air+Oxygen

scfh
lbm/min
kg/min
g/sec
g/min
mg/sec

Oxygen %

Save **Cancel** **Main Display**

Flow settings without Oxygen % option available

Settings

Medium

Units

Air
Oxygen
Air+Oxygen

scfh
lbm/min
kg/min
g/sec
g/min
mg/sec

Oxygen %

Save **Cancel** **Main Display**

Flow settings with Oxygen % available

Settings

Oxygen %

22 %



Save

Cancel

Oxygen % adjustment setting

Resetting Zero

The zero can be manually reset with flows displayed between -0.3 to 0.3 liters/min. This reset is active while the unit remains on. To remove the current zero setting, the unit can be deactivated and reactivated to the factory default.

Note that this functionality will only be available if Flow Rate is set as the "Main Display" measurement (see Main Display setting section above)

To Enable the resetting zero functionality:

1. Tap the gear icon in the top right of the default screen to enter the settings menu
2. Tap the **Set Zero** section
3. Tap the check box next to "Enable Re-Zeroing" (see **Fig. 3** below)
4. Return to the main screen to manually set the zero




Fig. 3
Default screen

To manually set the zero:

1. Tap the "Set Zero" icon at the top of the screen.
2. The flow rate numbers will turn red, indicating the zeroing functionality is activated.
3. Once the numbers are black, the flow rate should be zeroed.

Repeat the above steps until the flow rate is zero.

